

One Injection, **One Syringe**

PT Oneject Indonesia Sentul Factory:

Jl. Olympic Raya No.B9, Kawasan Industri Sentul, Bogor 16810, Indonesia

Bekasi Factory:

Jl. Anggrek VII No.41, Kawasan Industri Terpadu Indonesia China (KITIC), Kab. Bekasi 17337, Indonesia


Tel.: +62 21 8791 7422 (Hunting) Fax.: +62 21 8790 2761 E-mail: marketing@oneject.co.id





















Auto Disable Syringe

Every man, woman & child deserves the right to a safe injection

www.oneject.co.id

Oneject Auto Disable Syringe

ONEJECT Auto Disable Syringe design uses a simple yet effective technology by Star Syringe Limited (UK). After injection, the plunger will lock and immobilize automatically making the syringe unusable. This design is WHO / PQS approved and certified to all Auto Disable Syinges from PT Oneject Indonesia.





- 1. Auto Disable Syringe can only be used once.
- 2. After injection, plunger automatically locked and will break
- 3. Minimizes disease contamination from syringes.
- 4. Minimizes syringes misuse such as use more than once.
- 5. Protect Patients and Healthcare Workers from contamination.
- 6. Increase patient trust in Healthcare Workers.

Oneject AD Syringe 0.05ml 26G x 3/8" (10 mm) POS-WHO No: E008/030



Oneject AD Syringe 0.5ml 24G x 1" (25 mm) PQS-WHO No: E008/019



Oneject Hypodermic Syringe with Reuse Prevention Feature 2 ml 23G x 1 1/4" (32 mm)

PQS-WHO No:



Oneject Hypodermic Syringe with Reuse Prevention Feature 5 ml

21G x 1 1/2" (38 mm) POS-WHO No: E013/024



Oneject AD Syringe 0.3ml 24G x 1" (25 mm)



Oneject Hypodermic Syringe with Reuse Prevention Feature 1 ml

23G x 1" (25 mm) PQS-WHO No: E013/023



Oneject Hypodermic Syringe with Reuse Prevention Feature

23G x 1 1/4" (32 mm) POS-WHO No: E013/019



one injection, one syringe

